ADVANCE SCIENCE

INNOVATION | PRECISION | PERFORMANCE

SEE SAW SHAKER

MODEL NO. AS-SSS-01

- GENTLE ROCKING MOTION
- STABLE PERFORMANCE
- EASY OPERATION
- DIGITAL DISPLAY AND CONTROL
- COMPACT DESIGN



PRODUCT INTRODUCTION

Advance Science See Saw Shaker is designed for gentle and efficient mixing of liquids in culture dishes, staining trays and other laboratory containers. It provides a smooth and consistent tilting motion at a fixed angle, ensuring thorough mixing with minimal splashing.
With digital control, reliable performance and a compact design, it is suitable for various laboratory applications in research, clinical and diagnostic laboratories.

KEY FEATURES

- Gentle and uniform tilting motion
- Digital control with LCD display
- Fixed tilt angle for consistent mixing
- Suitable for various containers and trays
- Stable and quiet operation
- Compact and space-saving design
- Durable construction for long-term use

TECHNICAL SPECIFICATIONS

Parameter	Specification
Pan/Plate Material	Aluminum
Tilt Angle in	9 degree
Shaking Speed Range (Seesaw speed) in RPM	41 - 50
Timer Available	Yes
Maximum Load Capacity	10 kilogram
Display	LCD

APPLICATIONS

- Mixing of cell culture media
- Staining and de-staining applications
- Gentle mixing of samples
- Laboratory research
- Clinical and diagnostic laboratories
- General laboratory applications



FRONT VIEW



LEFT VIEW



RIGHT VIEW

info@advancescience.in | www.advancescience.in

ADVANCE SCIENCE

MSME ISO CE